

# Product Information



## Prion Protein Monoclonal Antibody - 12F10

Cat No: A03221 - 200 µg

### General Data

<b>Shipping:</b>	Dry Ice
<b>Formulation:</b>	lyophilized IgG with BSA
<b>Host:</b>	Mouse
<b>Antigen:</b>	This anti-Prion Protein Monoclonal Antibody was raised against a human recombinant protein.
<b>Clone:</b>	12F10
<b>Isotype:</b>	IgG2a



**Application(s):** Reconstitute the content of the vial in 1 mL of water.

For immunohistochemistry, the optimal working dilution must be determined empirically (currently between 0.1 and 1 µg/mL). For western blot analysis of PrPc, dilute the antibody to a final concentration of 1 µg/mL.

**Specificity:** PrPc  
(+/-) Hamster, (+/-) Mouse, (+) Bovine, (+) Ovine, (+) Human

### Product Overview

Prion Protein (PrP) and namely its abnormal isoform, partially resistant to proteinase K (PrPres), is the only specific molecular marker of the Transmissible Spongiform Encephalopathies (TSEs) such as Bovine Spongiform Encephalopathie (BSE) or its human form, the New Variant of Creutzfeld-Jakob disease.

This antibody recognises the protein sequence within amino acids 142-160 (human numbering).

FP/18/23

For research laboratory use only – Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

#### Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac  
10 bis avenue Ampère  
78180 Montigny le Bretonneux - France  
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>  
<https://www.bertin-bioreagent.com>



## Prion Protein Monoclonal Antibody - 12F10

Cat No: A03221 - 200 µg

### Scientific Literature

Krasemann S. et al. Generation of monoclonal antibodies against prion proteins with an unconventional nucleic acid-based immunization strategy. *Journal of Biotechnology* 73, 119-129 (1999).

Krasemann S. et al. Induction of antibodies against human prion proteins (PrP) by DNA-mediated immunization of PrP 0/0 mice. *Journal of Immunological Methods* 199, 109-118 (1996).

Krasemann S. et al. Generation of monoclonal antibody against Human prion proteins in PrP 0/0 Mice.

*Molecular Medicine*, 2, 725-731 (1996).

FP/18/23

**For research laboratory use only – Not for human diagnostic use.**

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Bioreagent does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Bioreagent be liable for incidental, consequential or punitive damages.

#### Contact Bertin Bioreagent

Parc d'Activités du Pas du Lac  
10 bis avenue Ampère  
78180 Montigny le Bretonneux - France  
Tel.: +33 (0)139 306 036

<https://www.bertin-bioreagent.com/pa206/contact-us>  
<https://www.bertin-bioreagent.com>